

## Recombinant Human ANXA13

Catalog No: #GP10168



Package Size: #GP10168-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ANXA13
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 56-316 amino acids of human ANXA13
Target Name	annexin A13
Other Names	ISA; ANX13
Accession No.	Swissprot:P27216Gene Accession:BC125158
Uniprot	P27216
GeneID	312;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. The specific function of this gene has not yet been determined; however, it is associated with the plasma membrane of undifferentiated, proliferating endothelial cells and differentiated villus enterocytes. Alternatively spliced transcript variants encoding different isoforms have been identified.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only