

Recombinant Human KRT18

Catalog No: #GP10170



Package Size: #GP10170-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human KRT18
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 130-430 amino acids of human KRT18
Target Name	keratin 18, type I
Other Names	K18; CYK18
Accession No.	Swissprot:P05783Gene Accession:BC000180
Uniprot	P05783
GeneID	3875;
Storage	-20~-80°C, pH 7.6 PBS

Background

KRT18 encodes the type I intermediate filament chain keratin 18. Keratin 18, together with its filament partner keratin 8, are perhaps the most commonly found members of the intermediate filament gene family. They are expressed in single layer epithelial tissues of the body. Mutations in this gene have been linked to cryptogenic cirrhosis. Two transcript variants encoding the same protein have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only