

Recombinant Human KRT17

Catalog No: #GP10172



Package Size: #GP10172-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human KRT17
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 102-346 amino acids of human KRT17
Target Name	keratin 17
Other Names	PC; K17; PC2; 39.1; CK-17; PCHC1
Accession No.	Swissprot:Q04695Gene Accession:BC000159
Uniprot	Q04695
GeneID	3872;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes the type I intermediate filament chain keratin 17, expressed in nail bed, hair follicle, sebaceous glands, and other epidermal appendages. Mutations in this gene lead to Jackson-Lawler type pachyonychia congenita and steatocystoma multiplex.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only