

Recombinant Human KRT15

Catalog No: #GP10173



Package Size: #GP10173-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human KRT15
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 300 amino acids of human Keratin, type I cytoskeletal 15
Target Name	Keratin, type I cytoskeletal 15
Other Names	K15, CK15, K1CO
Accession No.	Swissprot:P19012Gene Accession:BC002641
Uniprot	P19012
GeneID	3866;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a member of the keratin gene family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. Most of the type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains and are clustered in a region on chromosome 17q21.2.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only