

# Recombinant Human KRT4

Catalog No: #GP10174



Package Size: #GP10174-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human KRT4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 151-461 amino acids of human KRT4
Target Name	keratin 4, type II
Other Names	K4; CK4; CK-4; CYK4; WSN1
Accession No.	Swissprot:P19013Gene Accession:BC042174
Uniprot	P19013
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This type II cytokeratin is specifically expressed in differentiated layers of the mucosal and esophageal epithelia with family member KRT13.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only