

Recombinant Human KRT31

Catalog No: #GP10180



Package Size: #GP10180-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human KRT31
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 100-400 amino acids of human Keratin-31
Target Name	Keratin-31
Other Names	HA1, Ha-1, hHa1, KRTHA1
Accession No.	Swissprot:Q15323Gene Accession:BC114468
Uniprot	Q15323
GeneID	3881;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a member of the keratin gene family. As a type I hair keratin, it is an acidic protein which heterodimerizes with type II keratins to form hair and nails. The type I hair keratins are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only