

## YY1AP1 Antibody

Catalog No: #42866

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	YY1AP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total YY1AP1 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human YY1AP1
Target Name	YY1AP1
Other Names	YAP; HCCA1; HCCA2; YY1AP
Accession No.	Swiss-Prot#: Q9H869 Gene ID: 55249
Uniprot	Q9H869
GeneID	55249;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Colorectal tissue using #42866 at dilution 1/100,



Immunohistochemical analysis of paraffin-embedded Human Thyroid tissue using #42866 at dilution 1/100,

## Background

The encoded gene product presumably interacts with YY1 protein; however, its exact function is not known. Alternative splicing results in multiple transcript variants encoding different isoforms.?

Note: This product is for in vitro research use only