

Recombinant Human ARF6

Catalog No: #GP10190



Package Size: #GP10190-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ARF6
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ADP-ribosylation factor 6
Accession No.	Swissprot:P62330Gene Accession:BC008918
Uniprot	P62330
GeneID	382;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the human ARF gene family, which is part of the RAS superfamily. The ARF genes encode small guanine nucleotide-binding proteins that stimulate the ADP-ribosyltransferase activity of cholera toxin and play a role in vesicular trafficking and as activators of phospholipase D. The product of this gene is localized to the plasma membrane, and regulates vesicular trafficking, remodelling of membrane lipids, and signaling pathways that lead to actin remodeling. A pseudogene of this gene is located on chromosome 7.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only