

Recombinant Human RAP1A

Catalog No: #GP10194



Package Size: #GP10194-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RAP1A
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	RAP1A, member of RAS oncogene family
Other Names	RAP1; C21KG; G-22K; KREV1; KREV-1; SMGP21
Accession No.	Swissprot:P62834Gene Accession:BC014086
Uniprot	P62834
GeneID	5906;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the Ras family of small GTPases. The encoded protein undergoes a change in conformational state and activity, depending on whether it is bound to GTP or GDP. This protein is activated by several types of guanine nucleotide exchange factors (GEFs), and inactivated by two groups of GTPase-activating proteins (GAPs). The activation status of the encoded protein is therefore affected by the balance of intracellular levels of GEFs and GAPs.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only