

Recombinant Human RAP1B

Catalog No: #GP10195



Package Size: #GP10195-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RAP1B
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	Ras-related protein Rap-1b
Other Names	K-REV, RAL1B
Accession No.	Swissprot:P61224Gene Accession:BC000176
Uniprot	P61224
GeneID	5908;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the RAS-like small GTP-binding protein superfamily. Members of this family regulate multiple cellular processes including cell adhesion and growth and differentiation. This protein localizes to cellular membranes and has been shown to regulate integrin-mediated cell signaling. This protein also plays a role in regulating outside-in signaling in platelets. Alternate splicing results in multiple transcript variants. Pseudogenes of this gene are found on chromosomes 3, 5, 6 and 9.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only