

# Recombinant Human GST Tag

Catalog No: #GP10201



Package Size: #GP10201-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human GST Tag
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	Glutathione S-Transferase
Accession No.	Swissprot:Gene Accession:
Storage	-20~-80°C, pH 7.6 PBS

## Background

GST (Glutathione S-Transferase) is a 26kDa protein which was widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins. GST isoenzymes appear to play a central role in the parasite detoxification system. Other functions are also suspected including a role in increasing the solubility of haematin in the parasite gut.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only