

Recombinant Human BSA

Catalog No: #GP10202



Package Size: #GP10202-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BSA
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-607 amino acids of Bovine BSA
Target Name	Bovine serum albumin
Accession No.	Swissprot:P02769Gene Accession:NP_851335
Uniprot	P02769
GeneID	280717;
Storage	-20~-80°C, pH 7.6 PBS

Background

Serum albumin, the main protein of plasma, has a good binding capacity for water, Ca²⁺, Na⁺, K⁺, fatty acids, hormones, bilirubin and drugs. Its main function is the regulation of the colloidal osmotic pressure of blood. Major zinc transporter in plasma, typically binds about 80% of all plasma zinc.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only