

Recombinant Human ADAMTS15

Catalog No: #GP10203



Package Size: #GP10203-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ADAMTS15
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human ADAM metalloproteinase with thrombospondin type 1 motif, 15
Target Name	ADAM metalloproteinase with thrombospondin type 1 motif, 15
Accession No.	Swissprot:Q8TE58Gene Accession:BC109114
Uniprot	Q8TE58
GeneID	170689;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the ADAMTS (a disintegrin and metalloproteinase with thrombospondin motifs) protein family. ADAMTS family members share several distinct protein modules, including a propeptide region, a metalloproteinase domain, a disintegrin-like domain, and a thrombospondin type 1 (TS) motif. Individual members of this family differ in the number of C-terminal TS motifs, and some have unique C-terminal domains. The protein encoded by this gene has a high sequence similarity to the proteins encoded by ADAMTS1 and ADAMTS8.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only