

Recombinant Human BIRC5

Catalog No: #GP10220



Package Size: #GP10220-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BIRC5
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	Baculoviral IAP repeat-containing protein 5
Other Names	API4, EPR-1
Accession No.	Swissprot:O15392Gene Accession:BC008718
Uniprot	O15392
GeneID	332;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is a member of the inhibitor of apoptosis (IAP) gene family, which encode negative regulatory proteins that prevent apoptotic cell death. IAP family members usually contain multiple baculovirus IAP repeat (BIR) domains, but this gene encodes proteins with only a single BIR domain. The encoded proteins also lack a C-terminus RING finger domain. Gene expression is high during fetal development and in most tumors, yet low in adult tissues. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only