

# Recombinant Human BMP4

Catalog No: #GP10221



Package Size: #GP10221-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human BMP4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 116 amino acids of human BMP4
Target Name	bone morphogenetic protein 4
Other Names	ZYME; BMP2B; OFC11; BMP2B1; MCOPS6
Accession No.	Swissprot:P12644Gene Accession:BC020546
Uniprot	P12644
GeneID	652;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a member of the bone morphogenetic protein family which is part of the transforming growth factor-beta superfamily. The superfamily includes large families of growth and differentiation factors. Bone morphogenetic proteins were originally identified by an ability of demineralized bone extract to induce endochondral osteogenesis in vivo in an extraskelatal site.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only