

Recombinant Human BTG1

Catalog No: #GP10225



Package Size: #GP10225-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BTG1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	B-cell translocation gene 1, anti-proliferative
Accession No.	Swissprot:P62324Gene Accession:BC016759
Uniprot	P62324
GeneID	694;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is a member of an anti-proliferative gene family that regulates cell growth and differentiation. Expression of this gene is highest in the G0/G1 phases of the cell cycle and downregulated when cells progressed through G1. The encoded protein interacts with several nuclear receptors, and functions as a coactivator of cell differentiation. This locus has been shown to be involved in a t(8;12)(q24;q22) chromosomal translocation in a case of B-cell chronic lymphocytic leukemia.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only