

Recombinant Human BTG4

Catalog No: #GP10226



Package Size: #GP10226-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BTG4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 27-206 amino acids of human BTG4
Target Name	BTG family member 4
Other Names	PC3B
Accession No.	Swissprot:Q9NY30Gene Accession:BC031045
Uniprot	Q9NY30
GeneID	54766;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a member of the BTG/Tob family. This family has structurally related proteins that appear to have antiproliferative properties. This encoded protein can induce G1 arrest in the cell cycle.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only