

Recombinant Human BTG3

Catalog No: #GP10227



Package Size: #GP10227-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BTG3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 22-252 amino acids of human BTG3
Target Name	BTG anti-proliferation factor 3
Other Names	ANA; TOB5; TOFA; APRO4; TOB55; ANA/BTG3
Accession No.	Swissprot:Q14201Gene Accession:BC011957
Uniprot	Q14201
GeneID	10950;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a member of the BTG/Tob family. This family has structurally related proteins that appear to have antiproliferative properties. This encoded protein might play a role in neurogenesis in the central nervous system. Two transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only