

Recombinant Human BTG2

Catalog No: #GP10228



Package Size: #GP10228-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BTG2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	BTG anti-proliferation factor 2
Other Names	PC3; APRO1; TIS21
Accession No.	Swissprot:P78543Gene Accession:BC105949
Uniprot	P78543
GeneID	7832;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a member of the BTG/Tob family. This family has structurally related proteins that appear to have antiproliferative properties. This encoded protein is involved in the regulation of the G1/S transition of the cell cycle.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only