

ZFAND3 Antibody

Catalog No: #42870

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZFAND3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZFAND3 protein.
Immunogen Description	Full length fusion protein of human ZFAND3
Target Name	ZFAND3
Other Names	TEX27
Accession No.	Swiss-Prot#: Q9H8U3Gene ID: 60685
Uniprot	Q9H8U3
GeneID	60685;
Concentration	1.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Liver cancer tissue using #42870 at dilution 1/100.

Background

ZFAND3 contains 1 A20-type zinc finger and 1 AN1-type zinc finger. The exact function of ZFAND3 remains unknown.

Note: This product is for in vitro research use only