

# Recombinant Human CNN1

Catalog No: #GP10237



Package Size: #GP10237-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CNN1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	calponin 1
Other Names	SMCC; Sm-Calp; HEL-S-14
Accession No.	Swissprot:P51911Gene Accession:BC022015
Uniprot	P51911
GeneID	1264;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Thin filament-associated protein that is implicated in the regulation and modulation of smooth muscle contraction. It is capable of binding to actin, calmodulin, troponin C and tropomyosin. The interaction of calponin with actin inhibits the actomyosin Mg-ATPase activity (By similarity).

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only