

Recombinant Human CAMK1G

Catalog No: #GP10239



Package Size: #GP10239-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CAMK1G
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 23-277 amino acids of human CAMK1G
Target Name	calcium/calmodulin-dependent protein kinase IG
Other Names	VWS1; CLICK3; CLICKIII; dJ272L16.1
Accession No.	Swissprot:Q96NX5Gene Accession:BC032787
Uniprot	Q96NX5
GeneID	57172;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein similar to calcium/calmodulin dependent protein kinase, however, its exact function is not known.

Calcium/calmodulin-dependent protein kinase belonging to a proposed calcium-triggered signaling cascade. In vitro phosphorylates transcription factor CREB1 (By similarity).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only