

Recombinant Human CAP1

Catalog No: #GP10241



Package Size: #GP10241-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CAP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 319-453 amino acids of human Adenylyl cyclase-associated protein 1
Target Name	Adenylyl cyclase-associated protein 1
Other Names	CAP, CAP1-PEN
Accession No.	Swissprot:Q01518Gene Accession:BC013963
Uniprot	Q01518
GeneID	10487;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is related to the *S. cerevisiae* CAP protein, which is involved in the cyclic AMP pathway. The human protein is able to interact with other molecules of the same protein, as well as with CAP2 and actin. Alternatively spliced transcript variants have been identified.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only