

Recombinant Human CAP2

Catalog No: #GP10242



Package Size: #GP10242-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CAP2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a C terminal 300 amino acids of human CAP, adenylate cyclase-associated protein, 2 (yeast)
Target Name	CAP, adenylate cyclase-associated protein, 2 (yeast)
Accession No.	Swissprot:P40123Gene Accession:BC008481
Uniprot	P40123
GeneID	10486;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene was identified by its similarity to the gene for human adenylyl cyclase-associated protein. The function of the protein encoded by this gene is unknown. However, the protein appears to be able to interact with adenylyl cyclase-associated protein and actin.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only