

Recombinant Human CA10

Catalog No: #GP10249



Package Size: #GP10249-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CA10
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	carbonic anhydrase X
Other Names	CARPX; CA-RPX; HUCEP-15
Accession No.	Swissprot:Q9NS85Gene Accession:BC020577
Uniprot	Q9NS85
GeneID	56934;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that belongs to the carbonic anhydrase family of zinc metalloenzymes, which catalyze the reversible hydration of carbon dioxide in various biological processes. The protein encoded by this gene is an acatalytic member of the alpha-carbonic anhydrase subgroup, and it is thought to play a role in the central nervous system, especially in brain development. Multiple transcript variants encoding the same protein have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only