

Recombinant Human CBX2

Catalog No: #GP10263



Package Size: #GP10263-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CBX2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	chromobox homolog 2
Other Names	M33, CDCA6, SRXY5
Accession No.	Swissprot:Q14781Gene Accession:BC004252
Uniprot	Q14781
GeneID	84733;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a component of the polycomb multiprotein complex, which is required to maintain the transcriptionally repressive state of many genes throughout development via chromatin remodeling and modification of histones. Disruption of this gene in mice results in male-to-female gonadal sex reversal. Mutations in this gene are also associated with gonadal dysgenesis in humans. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only