

Recombinant Human CCK

Catalog No: #GP10264



Package Size: #GP10264-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CCK
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 21-115 amino acids of human CCK
Target Name	cholecystokinin
Accession No.	Swissprot:P06307Gene Accession:BC008283
Uniprot	P06307
GeneID	885;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the gastrin/cholecystokinin family of proteins. The encoded preproprotein is proteolytically processed to generate multiple protein products, including the peptide hormones cholecystokinin-8, -12, -33, and others. The encoded peptides have been shown to regulate gastric acid secretion and food intake. A sulfated form of cholecystokinin-8 may modulate neuronal activity in the brain. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only