

## ZFP2 Antibody

Catalog No: #42874

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZFP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZFP2 protein.
Immunogen Description	Fusion protein of human ZFP2
Target Name	ZFP2
Other Names	ZNF751
Accession No.	Swiss-Prot#: Q6ZN57 Gene ID: 80108
Uniprot	Q6ZN57
GeneID	80108;
Concentration	1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Breast cancer tissue using #42874 at dilution 1/100,

## Background

Probable transcription factor involved in neuronal differentiation and/or phenotypic maintenance.

---

Note: This product is for in vitro research use only