

# Recombinant Human CXCL10

Catalog No: #GP10272



Package Size: #GP10272-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CXCL10
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 22-98 amino acids of human CXCL10
Target Name	C-X-C motif chemokine ligand 10
Other Names	C7; IFI10; INP10; IP-10; crg-2; mob-1; SCYB10; gIP-10
Accession No.	Swissprot:P02778Gene Accession:BC010954
Uniprot	P02778
GeneID	3627;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This antimicrobial gene encodes a chemokine of the CXC subfamily and ligand for the receptor CXCR3. Binding of this protein to CXCR3 results in pleiotropic effects, including stimulation of monocytes, natural killer and T-cell migration, and modulation of adhesion molecule expression.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only