

Recombinant Human CDK20

Catalog No: #GP10273



Package Size: #GP10273-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CDK20
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	cyclin-dependent kinase 20
Other Names	P42, CCRK, CDCH, PNQALRE
Accession No.	Swissprot:Q8IZL9Gene Accession:BC002655
Uniprot	Q8IZL9
GeneID	23552;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene contains a kinase domain most closely related to the cyclin-dependent protein kinases. The encoded kinase may activate cyclin-dependent kinase 2 and is involved in cell growth. Alternatively spliced transcript variants encoding distinct isoforms have been reported.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only