

## Recombinant Human CDC20

Catalog No: #GP10274



Package Size: #GP10274-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CDC20
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 182-480 amino acids of human Cell division cycle protein 20 homolog
Target Name	Cell division cycle protein 20 homolog
Other Names	CDC20A, p55CDC, bA276H19.3
Accession No.	Swissprot:Q12834Gene Accession:BC010044
Uniprot	Q12834
GeneID	991;
Storage	-20~-80°C, pH 7.6 PBS

## Background

CDC20 appears to act as a regulatory protein interacting with several other proteins at multiple points in the cell cycle. It is required for two microtubule-dependent processes, nuclear movement prior to anaphase and chromosome separation.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only