

Recombinant Human CEACAM3

Catalog No: #GP10287

Package Size: #GP10287-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CEACAM3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 45-155 amino acids of human carcinoembryonic antigen-related cell adhesion molecule 3
Target Name	carcinoembryonic antigen-related cell adhesion molecule 3
Other Names	CEA; CGM1; W264; W282; CD66D
Accession No.	Swissprot:P40198Gene Accession:BC106728
Uniprot	P40198
GenID	1084;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the family of carcinoembryonic antigen-related cell adhesion molecules (CEACAMs), which are used by several bacterial pathogens to bind and invade host cells. The encoded transmembrane protein directs phagocytosis of several bacterial species that is dependent on the small GTPase Rac. It is thought to serve an important role in controlling human-specific pathogens by the innate immune system. Alternatively spliced transcript variants have been described.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only