

Recombinant Human CEACAM21

Catalog No: #GP10289



Package Size: #GP10289-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CEACAM21
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 35-191 amino acids of human carcinoembryonic antigen-related cell adhesion molecule 21
Target Name	carcinoembryonic antigen-related cell adhesion molecule 21
Other Names	CEACAM3, R29124_1
Accession No.	Swissprot:Q3KPI0Gene Accession:BC012001
Uniprot	Q3KPI0
GeneID	90273;
Storage	-20~-80°C, pH 7.6 PBS

Background

CEACAM21 is a gene which encodes carcinoembryonic antigen. It belongs to the immunoglobulin superfamily-CEA family. Three isoforms of the human protein are produced by alternative splicing. High levels of expression in hematopoietic system, brain and muscle.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only