

Recombinant Human CEBPE

Catalog No: #GP10291



Package Size: #GP10291-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CEBPE
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	CCAAT/enhancer binding protein (C/EBP), epsilon
Other Names	CRP1; C/EBP-epsilon
Accession No.	Swissprot:Q15744Gene Accession:BC035797
Uniprot	Q15744
GeneID	1053;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a bZIP transcription factor which can bind as a homodimer to certain DNA regulatory regions. It can also form heterodimers with the related protein CEBP-delta. The encoded protein may be essential for terminal differentiation and functional maturation of committed granulocyte progenitor cells. Mutations in this gene have been associated with Specific Granule Deficiency, a rare congenital disorder. Multiple variants of this gene have been described, but the full-length nature of only one has been determined.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only