

Recombinant Human CMTM2

Catalog No: #GP10299



Package Size: #GP10299-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CMTM2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 115 amino acids of human CKLF-like MARVEL transmembrane domain containing 2
Target Name	CKLF-like MARVEL transmembrane domain containing 2
Other Names	CKLFSF2
Accession No.	Swissprot:Q8TAZ6Gene Accession:BC025354
Uniprot	Q8TAZ6
GeneID	146225;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene belongs to the chemokine-like factor gene superfamily, a novel family that links the chemokine and the transmembrane 4 superfamilies of signaling molecules. The protein encoded by this gene may play an important role in testicular development.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only