

Recombinant Human CLEC4D

Catalog No: #GP10303



Package Size: #GP10303-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CLEC4D
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 39-215 amino acids of human CLEC4D
Target Name	C-type lectin domain family 4 member D
Other Names	MCL; MPCL; CD368; CLEC6; CLEC-6; CLECSF8; Dectin-3
Accession No.	Swissprot:Q8WXI8Gene Accession:BC032313
Uniprot	Q8WXI8
GeneID	338339;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the C-type lectin/C-type lectin-like domain (CTL/CTLD) superfamily. Members of this family share a common protein fold and have diverse functions, such as cell adhesion, cell-cell signalling, glycoprotein turnover, and roles in inflammation and immune response. This gene is closely linked to other CTL/CTLD superfamily members on chromosome 12p13 in the natural killer gene complex region.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only