

Recombinant Human CLEC2D

Catalog No: #GP10304



Package Size: #GP10304-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CLEC2D
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 60-155 amino acids of human C-type lectin domain family 2, member D
Target Name	C-type lectin domain family 2, member D
Other Names	CLAX; LLT1; OCIL
Accession No.	Swissprot:Q9UHP7Gene Accession:BC019883
Uniprot	Q9UHP7
GeneID	29121;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the natural killer cell receptor C-type lectin family. The encoded protein inhibits osteoclast formation and contains a transmembrane domain near the N-terminus as well as the C-type lectin-like extracellular domain. Several alternatively spliced transcript variants have been identified for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only