

RNF115 Antibody

Catalog No: #42878

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RNF115 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total RNF115 protein.
Immunogen Type	protein
Immunogen Description	Full length fusion protein of human RNF115
Target Name	RNF115
Other Names	BCA2; ZNF364
Accession No.	Swiss-Prot#: Q9Y4L5 Gene ID: 27246
Uniprot	Q9Y4L5
GeneID	27246;
Concentration	0.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Esophagus cancer tissue using #42878 at dilution 1/100.

Background

E3 ubiquitin-protein ligase that mediates E2-dependent, 'Lys-48'- and/or 'Lys-63'-linked polyubiquitination of substrates and may play a role in diverse biological processes. Through their polyubiquitination, may play a role in the endosomal trafficking and degradation of membrane receptors including EGFR, FLT3, MET and CXCR4.

Note: This product is for in vitro research use only