

Recombinant Human CNTFR

Catalog No: #GP10310



Package Size: #GP10310-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CNTFR
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 23-342 amino acids of human CNTFR
Target Name	ciliary neurotrophic factor receptor
Accession No.	Swissprot:P26992Gene Accession:BC001492
Uniprot	P26992
GeneID	1271;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the type 1 cytokine receptor family. The encoded protein is the ligand-specific component of a tripartite receptor for ciliary neurotrophic factor, which plays a critical role in neuronal cell survival, differentiation and gene expression. Binding of ciliary neurotrophic factor to the encoded protein recruits the transmembrane components of the receptor, gp130 and leukemia inhibitory factor receptor, facilitating signal transduction.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only