

Recombinant Human CDX1

Catalog No: #GP10320



Package Size: #GP10320-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CDX1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	caudal type homeobox 1
Accession No.	Swissprot:P47902Gene Accession:BC096252
Uniprot	P47902
GeneID	1044;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is a member of the caudal-related homeobox transcription factor gene family. The encoded DNA-binding protein regulates intestine-specific gene expression and enterocyte differentiation. It has been shown to induce expression of the intestinal alkaline phosphatase gene, and inhibit beta-catenin/T-cell factor transcriptional activity.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only