

# Recombinant Human CRMP1

Catalog No: #GP10322



Package Size: #GP10322-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CRMP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 249 amino acids of human collapsin response mediator protein 1
Target Name	collapsin response mediator protein 1
Other Names	DRP1, DRP-1, CRMP-1, DPYSL1, ULIP-3
Accession No.	Swissprot:Q14194Gene Accession:BC000252
Uniprot	Q14194
GeneID	1400;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of a family of cytosolic phosphoproteins expressed exclusively in the nervous system. The encoded protein is thought to be a part of the semaphorin signal transduction pathway implicated in semaphorin-induced growth cone collapse during neural development. Alternative splicing results in multiple transcript variants.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only