

Recombinant Human CRP

Catalog No: #GP10324



Package Size: #GP10324-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CRP
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-91 amino acids of human C-reactive protein, pentraxin-related
Target Name	C-reactive protein, pentraxin-related
Other Names	PTX1
Accession No.	Swissprot:P02741Gene Accession:BC020766
Uniprot	P02741
GenID	1401;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene belongs to the pentaxin family. It is involved in several host defense related functions based on its ability to recognize foreign pathogens and damaged cells of the host and to initiate their elimination by interacting with humoral and cellular effector systems in the blood. Consequently, the level of this protein in plasma increases greatly during acute phase response to tissue injury, infection, or other inflammatory stimuli.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only