

Recombinant Human GJA5

Catalog No: #GP10333



Package Size: #GP10333-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GJA5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 227-358 amino acids of human GJA5
Target Name	gap junction protein alpha 5
Other Names	CX40; ATFB11
Accession No.	Swissprot:P36382Gene Accession:BC013313
Uniprot	P36382
GeneID	2702;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is a member of the connexin gene family. The encoded protein is a component of gap junctions, which are composed of arrays of intercellular channels that provide a route for the diffusion of low molecular weight materials from cell to cell. Mutations in this gene may be associated with atrial fibrillation. Alternatively spliced transcript variants encoding the same isoform have been described.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only