

Recombinant Human CYCS

Catalog No: #GP10347



Package Size: #GP10347-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CYCS
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	cytochrome c, somatic
Other Names	CYC; HCS; THC4
Accession No.	Swissprot:P99999Gene Accession:BC015130
Uniprot	P99999
GeneID	54205;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a small heme protein that functions as a central component of the electron transport chain in mitochondria. The encoded protein associates with the inner membrane of the mitochondrion where it accepts electrons from cytochrome b and transfers them to the cytochrome oxidase complex. This protein is also involved in initiation of apoptosis. Mutations in this gene are associated with autosomal dominant nonsyndromic thrombocytopenia. Numerous processed pseudogenes of this gene are found throughout the human genome.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only