

Recombinant Human LAMP2

Catalog No: #GP10352



Package Size: #GP10352-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Recombinant Human LAMP2 |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 29-375 amino acids of human lysosomal-associated membrane protein 2 |
| Target Name | lysosomal-associated membrane protein 2 |
| Other Names | LAMPB; CD107b; LAMP-2; LGP110 |
| Accession No. | Swissprot:P13473Gene Accession:BC002965 |
| Uniprot | P13473 |
| GenID | 3920; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may play a role in tumor cell metastasis. It may also function in the protection, maintenance, and adhesion of the lysosome. Alternative splicing of this gene results in multiple transcript variants encoding distinct proteins

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only