

Recombinant Human DAP3

Catalog No: #GP10357

Package Size: #GP10357-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DAP3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a C terminal 300 amino acids of human death associated protein 3
Target Name	death associated protein 3
Other Names	DAP-3, MRPS29, MRP-S29, bMRP-10
Accession No.	Swissprot:P51398Gene Accession:BC107488
Uniprot	P51398
GeneID	7818;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a 28S subunit protein that also participates in apoptotic pathways which are initiated by tumor necrosis factor-alpha, Fas ligand, and gamma interferon. This protein potentially binds ATP/GTP and might be a functional partner of the mitoribosomal protein S27. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. Pseudogenes corresponding to this gene are found on chromosomes 1q and 2q.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only