

Recombinant Human DTX3

Catalog No: #GP10361



Package Size: #GP10361-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DTX3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 248 amino acids of human deltex homolog 3 (Drosophila)
Target Name	deltex homolog 3 (Drosophila)
Other Names	RNF154
Accession No.	Swissprot:Q8N9I9Gene Accession:BC114498
Uniprot	Q8N9I9
GeneID	196403;
Storage	-20~-80°C, pH 7.6 PBS

Background

DTX3 functions as an E3 ubiquitin ligase. It is regulator of Notch signaling, a signaling pathway involved in cell-cell communications that regulates a broad spectrum of cell-fate determinations. Probably acts both as a positive and negative regulator of Notch, depending on the developmental and cell context. Functions as an ubiquitin ligase protein in vitro, suggesting that it may regulate the Notch pathway via some ubiquitin ligase activity.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only