

Recombinant Human DKK3

Catalog No: #GP10364



Package Size: #GP10364-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DKK3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 147-350 amino acids of human dickkopf WNT signaling pathway inhibitor 3
Target Name	dickkopf WNT signaling pathway inhibitor 3
Other Names	RIG; REIC
Accession No.	Swissprot:Q9UBP4Gene Accession:BC007660
Uniprot	Q9UBP4
GenelD	27122;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that is a member of the dickkopf family. The secreted protein contains two cysteine rich regions and is involved in embryonic development through its interactions with the Wnt signaling pathway. The expression of this gene is decreased in a variety of cancer cell lines and it may function as a tumor suppressor gene. Alternative splicing results in multiple transcript variants encoding the same protein.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only