

# Recombinant Human DKKL1

Catalog No: #GP10365



Package Size: #GP10365-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human DKKL1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 31-242 amino acids of human DKKL1
Target Name	dickkopf like acrosomal protein 1
Other Names	SGY; CT34; SGY1; SGY-1
Accession No.	Swissprot:Q9QZL9Gene Accession:BC030581
Uniprot	Q9QZL9
GeneID	50722;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The dickkopf protein family interacts with the Wnt signaling pathway and its members are characterized by two conserved cysteine-rich domains. This gene encodes a secreted protein that has low sequence similarity to the dickkopf-3 protein. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only