

Recombinant Human DKK4

Catalog No: #GP10366



Package Size: #GP10366-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DKK4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 21-224 amino acids of human dickkopf WNT signaling pathway inhibitor 4
Target Name	dickkopf WNT signaling pathway inhibitor 4
Other Names	DKK-4
Accession No.	Swissprot:Q9UBT3Gene Accession:BC107047
Uniprot	Q9UBT3
GenelD	27121;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that is a member of the dickkopf family. The secreted protein contains two cysteine rich regions and is involved in embryonic development through its interactions with the Wnt signaling pathway. Activity of this protein is modulated by binding to the Wnt co-receptor and the co-factor kremen 2.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only